



Preva Radiation Dose Table

10/30/2006

Instruments Used

Radcal Model 9010	s/n	90-2215
Radcal Model 4083	s/n	43-0221
Ion Chamber Model 10	s/n	10294

Preva Tube Head	s/n	TA60046
-----------------	-----	---------

59mm round cone

measured kVp: measured with Radcal 4083
measured mA: read from Preva display
ion chamber placed at the end of the cone

All exposures are with 0.100 second exposure times.

kVp	mA	mR
60	4	34.5
60	5	42.6
60	6	48.6
60	7	56.0
65	4	44.1
65	5	53.3
65	6	61.5
65	7	68.1
70	4	51.8
70	5	61.5
70	6	70.0